

METAFECTENE® SI⁺ Technical Note

Dr. Natalie Falkenberg, Michaela Huber, Rebecca Mall, Prof. Dr. Michaela Aubele

Institute of Pathology, Hemholtz Zentrum München, Ingolstaedter Landstrasse 1,
85764 Neuherberg, Germany

Materials: ON-TARGET plus siRNA targeting human membrane-associated protein
(target protein 1)
ON-TARGET plus siRNA targeting human target protein 2
ON-TARGET plus siRNA targeting human GAPDH (positive control)
ON-TARGET plus non-targeting siRNA (negative control)
(all siRNAs from Thermo Scientific/ABgene, Epsom Surrey, UK and at a
concentration of 20 pmol/μl)
Two sterile 6-well culture plates
Sterile Eppendorf tubes
0.05 % Trypsin-EDTA solution
Complete medium: RPMI with Glutamax, 10 % fetal calf serum, 0.25 %
each of Penicillin/Streptomycin (all from Life Technologies, Darmstadt,
DE), 0,4 % bovine insulin (Sigma, St. Louis, MO, USA)
BT549 cells (Breast cancer)

Transfection reagents:

METAFECTENE® SI⁺

Lipofectamine® 2000 (Life Technologies, Darmstadt, DE)

OptiMEM (Life Technologies, Darmstadt, DE)

Antibiotic-free medium: RPMI with Glutamax, 10 % fetal calf serum (both from
Life Technologies, Darmstadt, DE), 0,4 % bovine insulin (Sigma, St. Louis,
MO, USA)

All steps were conducted under sterile conditions.

Protocol for the transient transfection of BT549 cells

A total of 0.5×10^6 cells/well was seeded and grown to 80 % confluence in complete medium at 37°C in a humidified incubator under 5 % CO₂ conditions. By the day of transfection, the complete medium was replaced by antibiotic-free medium and the cells were incubated for a minimum of 1 h.

In the meantime, the respective transfection reagents were pipetted into an 1,5 ml Eppendorf tube (Table 1 and 2).

For transfections using METAFECTENE® SI⁺:

All steps were performed according to manufacturer's instructions. 150 μl 1 x SI⁺ Buffer was mixed with METAFECTENE® SI⁺ and 200 pmol siRNA were added and carefully mixed (Table 1). After 15 min incubation at room temperature, the METAFECTENE® SI⁺ - siRNA complexes were pipetted onto the cells and mixed by swinging the plates for a few times.

	SI ⁺ Buffer	METAFECTENE [®] SI ⁺	siRNA
1. mock:	150 µl	7,2 µl	-
2. + siRNA.NK:	150 µl	7,2 µl	10 µl
3. + siRNA.GAPDH:	150 µl	7,2 µl	10 µl
4. + siRNA.TP1:	150 µl	7,2 µl	10 µl
5. + siRNA.TP2:	150 µl	7,2 µl	10 µl
6. + siRNA.TP1 + siRNA.TP2:	140 µl	14,4 µl	10 µl + 10 µl

Table 1: The amounts of respective reagents for the transfection using METAFECTENE[®] SI⁺. TP: target protein, NK: negative control.

For transfections using Lipofectamine[®] 2000:

All steps were performed according to manufacturer's instructions. The respective siRNAs were added to OptiMEM and in a separate tube, 10 µl of Lipofectamine[®] 2000 were mixed with OptiMEM and incubated for 5 min at room temperature (Table 2). Both approaches were combined, mixed carefully and incubated for 20 min at room temperature. Then, the Lipofectamine[®] 2000-siRNA complexes were pipetted onto the cells and mixed by swinging the plates for a few times.

	I.		II.	
	siRNA	OptiMEM	Lipofectamine [®] 2000	OptiMEM
1. mock:	-	250 µl	10 µl	240 µl
2. +siRNA.NK:	10 µl	240 µl	10 µl	240 µl
3. +siRNA.GAPDH:	10 µl	240 µl	10 µl	240 µl
4. +siRNA.TP1:	10 µl	240 µl	10 µl	240 µl
5. +siRNA.TP2:	10 µl	240 µl	10 µl	240 µl
6. +siRNA.TP1 +siRNA.TP2:	10 µl + 10 µl	230 µl	10 µl + 10 µl	230 µl

Table 2: The amounts of respective reagents for the transfection using Lipofectamine[®] 2000. TP: target protein, NK: negative control.

In both approaches, the cells were incubated for 6 h at 37°C in a humidified incubator under 5 % CO₂ conditions, then the supernatant was removed, the cells were washed with PBS, transferred into cell culture flasks and incubated in complete medium. 72 h post-transfection, the supernatant was removed, the cells were trypsinized, pelleted (1500 rpm, 5 min) and lysed in T-PER[®] tissue protein extraction reagent buffer (Thermo Scientific, Rockford, IL, USA) supplemented with Complete Mini Protease inhibitor (Roche, Mannheim, DE) and Phosphatase Inhibitor cocktail 2 and 3 (Sigma, St. Louis, MO, USA).

SDS-PAGE and Western blots:

The protein concentrations were determined using the BCA Protein Assay Reagent (Thermo Scientific, Rockford, IL, USA). For Western blot analysis, 50 µg of protein lysates were separated by SDS-PAGE and transferred onto a 45 µm nitrocellulose membrane (GE Healthcare/Amersham, Buckinghamshire, UK). The membrane was blocked in a 5 % skim milk T-BST solution for 1 h at 37°C. The antibodies against the target proteins and the controls were incubated overnight in 5% BSA/T-BST at 4°C and were detected using HRP-linked secondary antibodies incubating for 1 h in 5% skim milk solution at room temperature.




Figure 1: Western blot analysis of BT549 cell lysates derived from transient (72h) transfections using Lipofectamine® 2000 and METAFECTENE® SI⁺. TP: target protein, NK: negative control, GAPDH: positive control.

Conclusions

The Western blot analysis (Figure 1) demonstrates a very efficient downregulation of the target proteins using either METAFECTENE® SI⁺ or Lipofectamine® 2000. Although both approaches were conducted under the same conditions and at the same time, GAPDH (positive control) seems to be stronger downregulated when METAFECTENE® SI⁺ was used. Furthermore, in comparison with the Lipofectamine® 2000 protocol, the transfection protocol using METAFECTENE® SI⁺ is time-saving due to less pipetting and intermediate steps. This is very advantageous, in particular when a lot of approaches should be conducted. Moreover, no toxic effects on BT549 cells were observed using the transfection reagents under the previously described conditions.