

## **Transient transfection of GD25- $\beta_1$ mouse embryonal fibroblasts using METAFACTENE transfection reagent**

**C. Smerling, R.W. Farndale**

Department of Biochemistry, University of Cambridge, CB2 1QW, Cambridge,  
UK

### **Cell line**

GD25- $\beta_1$  cells are derived from the knock-out cell line GD25, which is deficient in integrin  $\beta_1$  expression<sup>1</sup>. GD25- $\beta_1$  cells have been stably re-transfected with the integrin  $\beta_1$  gene and are used as control cells.

### **Protocol**

The medium of cells at 90-100% apparent confluency was exchanged for DMEM with 10% FCS without antibiotics. For transient transfections in a 3.5 cm culture dish, 5  $\mu$ g of plasmid DNA and 10  $\mu$ l metafectene were mixed in a final volume of 200  $\mu$ l DMEM without FCS and antibiotics. This mixture was incubated for 15 – 20 min at room temperature, before being added to the cells. The cells were then incubated for 48 h with the DNA/lipid complexes before being subjected to flow cytometry.


For FACS analysis  $8 \times 10^5$  cells were detached with PBS / EDTA and resuspended in 400  $\mu$ l of PBS. 20  $\mu$ l of FITC-conjugated anti- $\alpha_2$  antibody were added and the cells were incubated in the dark on ice for 1h. The cells were washed three times with 1 ml PBS and then resuspended in 400  $\mu$ l of PBS ( $2 \times 10^6$  cells per ml). FACS analysis was done using a BD FACScan analyser at a wavelength of 515 – 545 nm.

### **Discussion**

We have established that using 5  $\mu$ g of plasmid DNA and 10  $\mu$ l metafectene per 3.5 cm dish results in a transfection efficiency of about 50% as determined by flow cytometry (Fig. 1a). This applies to the GFP-tagged aminoglycoside phosphotransferase, which is expressed from a SV40 promotor. The transmembrane protein integrin  $\alpha_2$  appears to be expressed in a percentage of cells slightly smaller than 50% (Fig. 1b). Both proteins are expressed from a plasmid derived from pCMV-Script<sup>®</sup> (Stratagene).

Next, we will attempt to establish a stable cell line with the plasmid used here for transient transfection and use METAFECTENE for screening mutants transiently expressed in these cells.

## Results



**Fig. 1:** Histograms of FACS analysis of transiently transfected GD25 cells 48 hours after transfection.

1. Fassler, R., *et al.* Lack of beta 1 integrin gene in embryonic stem cells affects morphology, adhesion, and migration but not integration into the inner cell mass of blastocysts. *J Cell Biol*, **128**(5): 979-88.