

Use of Metafectene for the Transfection of Mouse Embryonic Fibroblasts

Jihae Shin, Graduate Assistant, Dept. of Molecular Genetics and Microbiology,
University of Florida, FL USA

1. Transfection Protocol

MEFs were plated 24 hours prior to transfection and cultured in growth medium without antibiotics (DMEM + 15% FBS) at 37°C. For transfections in 6 well plates, 4µl/8µl /12µl of Metafectene was diluted into a total of 100µl of OptiMEM (Gibco) and 2µg of DNA was diluted into a total of 100µl of OptiMEM. The dilutions were mixed and incubated at room temperature for 30 minutes. Cells were washed with PBS and 1 ml of OptiMEM was added to each plate. Then Metafectene/DNA complex was added to plate dropwise. The transfections were incubated at 37. After 6 hours, the transfection mixture was aspirated and 2 ml of growth medium was added to each well. GFP signal was checked the next day.

2. Optimization of Metafectene Transfection

Primary mouse embryonic fibroblasts (MEF) were transfected with Metafectene under 3 different conditions in order to optimize the transfection efficiency. The ratio of DNA to Metafectene was varied. The condition of transfection is described below.

6-well plate (per well)

Cell type	Primary MEFs (P3)
Confluency	80-90%
DNA (pEGFP-C1)	2 µg
Metafectene	4µl / 8µl / 12µl

The transfections using 8µl or 12µl of Metafectene had **~40% GFP positive cells**.

3. Comparison of Transfection reagents

MEFs were also transfected with other lipid-based transfection reagents (e.g. Lipid L2). The same protocol was applied except the amount of Lipid L2 used.

6-well plate (per well)

Cell type	Primary MEFs (P3)
Confluency	80-90%
DNA (pEGFP-C1)	2 µg
Lipid L2	2µl / 4µl / 8µl

The transfections using Lipid L2 **had ~20% GFP positive cells**.

4. Conclusion

Metafectene can be used to transiently transfect primary MEFs. The reagent performed better than other lipid based transfection reagents tested. Usage of 4:1 or 6:1 Metafectene:DNA is recommended.