

## LRRK2-induced Neuronal Toxicity

Yulan Xiong/ Valina Dawson  
Neuroregeneration and Stem Cell Programs  
Department of Neurology  
Institute for Cell Engineering,  
Johns Hopkins University, School of Medicine  
733 N. Broadway, BRB 732  
Baltimore, Maryland 21205

### **Introduction:**

Mutations in the leucine-rich repeat kinase 2 (LRRK2) gene have been identified as a unambiguous cause of rare autosomal dominant forms of PD. However, the pathogenic role and associated biochemical pathways responsible for LRRK2-linked disease still remains unclear. To investigate the neurotoxicity induced by LRRK2, we transfected LRRK2 WT and G2019S mutant into mouse cortical neuronal culture by METAFECTENE EASY reagents.

### **Materials and methods:**

Transfect reagent: Metafectene EASY and Metafectene SI (Biontex, Munich, Germany);  
Plasmids: pcDNA3.1/MYC-HIS-LRRK2 –WT/G2019S, pcDNA3.1-eGFP;  
Cells: mouse cortical neuronal culture

### **Experimental procedures / transfection protocol:**

- 1). The mouse cortical neuronal culture grew on 24-well plates for one week.
- 2). Dilute 10XEASY buffer to 1X
- 3). Add 5ul METAFECTENE EASY in 100ul 1XEASY buffer and mix gently
- 4). Add 2.2ug DNA (2ug LRRK2, 0.2ug eGFP) into above mix, incubate for 15 mins at RT
- 5). Change the culture media to OPTI-MEM media (400ul per 24-well)
- 6). Add the mix from step 4 to the culture
- 7). Chang the media to neuronal growth media next morning
- 8). The viability of eGFP positive neurons using TUNEL staining or nuclear condensation/fragmentation by fluorescence microscopy were quantified to assess the neuronal toxicity.

### **Results and discussion:**


Transfection of 7-days old cortical neurons using Metafectene EASY resulted in a transfection efficiency around 1-2%, which is similar to the efficiency of Lipofectamine 2000 (Invitrogen, data not shown)

### **Conclusion / summary:**

Transfection efficiency was as high as compared to results achieved using transfection reagents of other manufacturers. We didn't see obvious decreased cytotoxicity compared to other transfection reagents.

Appendix: Tables and/or figures:

| Cell code | Primary | Class     | Species | Organ  | Reagent          | Growth Properties | Genetic Material | Efficiency | Toxicity |
|-----------|---------|-----------|---------|--------|------------------|-------------------|------------------|------------|----------|
| Neurons   | yes     | Mam-malia | mouse   | cortex | METAFECTENE EASY | adherent          | Plasmid          | 1-2%       | low      |



Empty vector+eGFP                      LRRK2-WT+eGFP                      LRRK2-G2019S+eGFP

**Figure 1. LRRK2 induced neuronal toxicity.** 7-days old mouse cortical neurons were co-transfected with LRRK2 and eGFP (10:1 ratio) by Metafectene Easy reagent. Viability was analyzed at 48 hrs post-transfection (DIV 9) with non-viable neurons exhibiting obvious neurite process and/or nuclear fragmentation. Representative images are showing eGFP postive neurons.