

## **Transfection rates of HEK 293 and NG 108-15 cells using METAFACTENE and METAFACTENE PRO**

Prof. Dr. Rüdiger Schulz and Andrea Wehmeyer  
Institute for Pharmacology, Toxicology und Pharmacy  
School of Veterinary Medicine, LMU Munich  
Königinstraße 16  
D - 80539 Munich, Germany

### **Experimental procedure:**

METAFACTENE and METAFACTENE PRO (Biontex Laboratories GmbH, Martinsried, Germany) were examined for transfection efficiency in HEK 293 and NG 108-15 permanent cell lines. The cells were seeded in 6-well-plates. The transfected material was DNA (in pcDNA3.1(+)) coding for the fusion protein  $\mu$ -opioid-receptor-EGFP (EGFP, enhanced green fluorescence protein). The transfection efficiency was judged by visually counting the number of fluorescent cells, neglecting the intensity of fluorescence. The experimental protocol considered various amounts of DNA and reagents, respectively. Data reflect the outcome of a single experimental set up conducted in duplicate.

DNA ( $\mu$ g)	Metafectene ( $\mu$ l)	HEK 293 %	NG 108-15 %	Metafectene Pro ( $\mu$ l)	HEK 293 %	NG 108-15 %
0.4	2	0	0	2	0	0
0.4	3	0	0	3	0	0
1.0	3	40	5-10	3	40	40
1.0	6	45	5-10	6	< 50	50
2.0	6	60	5-10	6	60	60
2.0	12	55	5-10	12	50	++
3.0	12	++	++	12	++	++
4.0	24	+++	+++	24	+++	+++

++: cells died

+++ : many cells died

### **Conclusion:**

The transfection efficiencies of METAFACTENE and METAFACTENE PRO are approximately identical concerning HEK 293 cells.

The transfection efficiency of METAFACTENE PRO is explicitly higher in NG 108-15 cells as compared to METAFACTENE.