

Transfection of *rab6* gene from *Leishmania donovani* using

Biont~~ex~~ K2 transfection system

Indira Singh Chauhan and Neeloo Singh*

Biochemistry Division, CSIR-Central Drug Research Institute, Jankipuram Extension, Sitapur

Road, Lucknow 226031, India

*Corresponding author. Tel: +91-9415002065; Fax: +91-522-22623405;

E-mail: neeloo888@yahoo.com

Materials and Methods

Cell culture


Promastigotes of *Leishmania donovani* strain (MHOM/IN/80/DD8) were maintained in Medium 199 (pH 7.2) (Sigma, St Louis, MO, USA) 2.05 mM L-glutamine, 12 mM HEPES buffer (Sigma), 10% (v/v) HIFBS, 100 U/mL penicillin and 100 µg/mL streptomycin. Promastigotes were grown in vented T25 tissue culture flasks and maintained at 25°C. Promastigote cultures were initiated at 10⁶ parasites/mL and subcultured every 3–4 day.

Cell transfection:

Late logarithmic phase of *Leishmania donovani* promastigotes (1x10⁷ cell/ml) were harvested by centrifugation at 4500 rpm for 5 min at room temperature. Promastigotes were treated with K2 multiplier and incubated for 2 hr at 25°C before transfection. For this K2 Multiplier (5 µL) was added slowly onto the culture media (250 µL) and mixed by gently. Transfection solution (lipoplex) was prepared in serum-free cell culture medium. 20 µg of plasmid DNA (pET28a+ and *Ldrab6* gene) was added to a 9 µL of K2 transfection solution and incubated for 20 min at room temperature. Transfection solution (lipoplex) was added to the promastigotes cells and incubated for 24 hr at 25°C. On the genomic level, transfection efficiency was monitored by amplified PCR product of *Ldrab6* gene encoded by the transfected plasmids.

Result:

39 *L. donovani* promastigotes were transfected with plasmid DNA (pET28a+ and *Ldrab6* gene)
40 using Biontex K2 transfection system. Transfection was confirmed by PCR using T7 promotor of
41 pET28a⁺ vector and *rab6* gene primer. Wild type promastigotes have 825 bp of *rab6* gene and
42 transfected promastigotes have 1020 bp of expressed *rab6* gene which was subcloned in
43 pET28a+ expression vector.



44

45 **Transfection of *Ldrab6* gene in *Leishmania donovani* by using K2 transfection system**
46 **(Biontex). After transfection, pellet was used for DNA isolation by Qiagen kit (Cat.**
47 **No.51104). Lane 1; PCR product of *rab6* gene from *L. donovani* (using *rab6* forward and**
48 **reverse primer) was used as control, Lane 2; PCR product of transfected *L. donovani* with**
49 ***Ldrab6* gene (using T7 promoter, Lot# 38530153 as forward primer and *rab6* reverse as**
50 **reverse primer).**

51

52 **Conclusions:**

53 K2 transfection system was used by us for the first time in *Leishmania donovani*. It is based on
54 cationic liposome formulation which is used to increase transfection efficiency of plasmid DNA.
55 Transfected parasites show a band of the expected size (~1020 bp) when compare with wild type
56 parasites (825 bp). The transfection rate was increased approximately 80% when compare with
57 other transfection reagents.